



AT TORBI 709M
AT TORBI 709MP preloaded
Technical Specifications



AT TORBI 709M

AT TORBI 709MP preloaded



AT TORBI® 709M

Optic Design	Monofocal, bitoric, aspheric (aberration neutral)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, bitoric, MICS
Incision Size	1.8 mm for IOLs from -10.0 to +28.0 D sphere with +1.0 to +4.0 D cylinder / 2.2 mm for all remaining diopter-cylinder-combinations up to +32.0 D
Company Labeled A-Constant ¹	118.3
Diopter Range	Sphere: -10.0 to +32.0 D, 0.5 D increments Cylinder: +1.0 to +12.0 D, 0.5 D increments
ACD	5.07
Implantation in	Capsular bag
Injector / Cartridge Set	For IOLs from -10.0 to +28.0 D sphere with +1.0 to +4.0 D cylinder*: AT.Shooter A2-2000 injector / AT.Smart Cartridge Set <i>or</i> AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set * For all remaining diopter-cylinder-combinations use VISCOJECT-BIO 2.2 Injector Set <i>or</i> AT.Shooter A1-2000 injector with VISCOJECT-BIO 2.2 Cartridge Set.

AT TORBI® 709MP preloaded

Optic Design	Monofocal, bitoric, aspheric (aberration neutral)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, bitoric, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	118.3
Diopter Range	Sphere: -10.0 to +28.0 D, 0.5 D increments Cylinder: +1.0 to +4.0 D, 0.5 D increments
ACD	5.07
Implantation in	Capsular bag
Injector / Cartridge Set ²	BLUEMIXS® 180

¹ Please refer to our web pages for optimized A-Constants. ² Please refer to our web pages for the most up-to-date references.



AT TORBI 709M / 709MP



BLUEMIXS 180 injector



AT.Shooter A1-2000
AT.Shooter A2-2000



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts



Carl Zeiss Meditec SAS
27 Avenue Paul Langevin
17180 Périgny - France
Mailing address:
Avenue Paul Langevin, BP5
17053 La Rochelle Cedex 9
France



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts

EN_32_022_0016111 Printed in Germany. CZ-11/2018 International edition: Only for sale in selected countries.
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representative for more information. Subject to changes in design and scope of delivery and as a result of ongoing technical development. AT TORBI and BLUEMIXS are registered trademarks of Carl Zeiss Meditec AG in Germany and/or other countries. VISCOJECT is a trademark of Meditel (CE0482).
© Carl Zeiss Meditec AG, 2018. All rights reserved.